

Soothe-Guard[®]

PROTECTIVE RADIOGRAPHY APRONS

Soothe-Guard AIR[™]

LEAD-FREE RADIOGRAPHY APRONS

SELECTION GUIDE



DENTSPLY
RINN

Soothe-Guard®

PROTECTIVE RADIOGRAPHY APRONS

+++++ – *The Dental Advisor*

The most comfortable aprons we've ever offered!

Unique, silky-smooth material has an exquisite feel that is both relaxing and comforting.

New collar design fits better for less gapping, greater patient comfort.

Sturdy **sewn-in hanging loops** for easy storage.

Soothe-Guard material front and back makes these aprons **easy to clean** (just wipe with damp cloth or disinfectant wipe), with **smooth draping**, and **exceptional durability**.



	Buttercup Yellow	Navy Blue	Dove Grey	Light Blue	Tan/Beige	Dusty Mauve
.3MM EQUIVALENCY						
Adult Apron with Collar	668050	669050	6610050	6611050	6612050	6613050
Adult Apron	668049	669049	6610049	6611049	6612049	6613049
Child Apron with Collar	668047	669047	6610047	6611047	6612047	6613047
Child Apron	668048	669048	6610048	6611048	6612048	6613048
Protective Collar	668051	669051	6610051	6611051	6612051	6613051
Pano-Cape	668053	669053	6610053	6611053	6612053	6613053

.5MM EQUIVALENCY						
Adult Apron with Collar	568050	569050	5610050	5611050	5612050	5613050
Adult Apron	568049	569049	5610049	5611049	5612049	5613049
Child Apron with Collar	568047	569047	5610047	5611047	5612047	5613047
Child Apron	568048	569048	5610048	5611048	5612048	5613048
Protective Collar	568051	569051	5610051	5611051	5612051	5613051
Pano-Cape	568053	569053	5610053	5611053	5612053	5613053



Buttercup Yellow

Navy Blue

Dove Grey

New!

Soothe-Guard AIR™

LEAD-FREE RADIOGRAPHY APRONS



All the features of our regular Soothe-Guard aprons, plus...

Innovative lead-free, alloy lining provides protection equivalent to lead with less weight. Go to Rinncorp.com for more info.

Up to **37% lighter** than lead-lined aprons.

Easier handling for dental professionals.

Environmentally friendly and requires **no special disposal procedure**.

	Buttercup Yellow	Navy Blue	Dove Grey	Light Blue	Tan/Beige	Dusty Mauve
.3MM EQUIVALENCY						
Adult Apron with Collar	868050	869050	8610050	8611050	8612050	8613050
Adult Apron	868049	869049	8610049	8611049	8612049	8613049
Child Apron with Collar	868047	869047	8610047	8611047	8612047	8613047
Child Apron	868048	869048	8610048	8611048	8612048	8613048
Protective Collar with Velcro	868051 868451	869051 869451	8610051	8611051	8612051	8613051
Pano-Cape	868053	869053	8610053	8611053	8612053	8613053
Adult Pano-Dual	868054	869054	8610054	8611054	8612054	8613054
Child Pano-Dual	868055	869055	8610055	8611055	8612055	8613055
Technician (Coat Style)	868056	869056	8610056	8611056	8612056	8613056
.5MM EQUIVALENCY						
Adult Apron with Collar	878050	879050				
Child Apron with Collar	878047	879047				
Adult Pano-Dual	878054	879054				



Light Blue

Tan/Beige

Dusty Mauve

Apron Styles



Not shown: Child PanoDual

Apron Dimensions and Weights

Description	Size		Weight							
	inches	millimeters	SOOTHE-GUARD AIR				SOOTHE-GUARD			
			.3mm Lead Equivalency		.5mm Lead Equivalency		.3mm Lead Equivalency		.5mm Lead Equivalency	
			lbs.	kgs.	lbs.	kgs.	lbs.	kgs.	lbs.	kgs.
Adult Patient	23.4 x 25.4	594 x 645	3.80	1.7	—	—	5.75	2.6	8.88	4.0
Adult Patient with Collar	23.0 x 25.0	584 x 635	4.05	1.8	6.60	2.99	6.50	3.0	9.56	4.4
Child Patient	20.9 x 19.7	531 x 500	2.80	1.3	—	—	3.81	1.7	5.88	2.7
Child Patient with Collar	20.5 x 19.6	521 x 498	2.70	1.2	4.00	1.81	4.50	2.0	7.06	3.2
Adult PanoDual	23.3 x 26.8	591 x 679	6.75	3.1	10.35	4.69	10.38	4.7	15.12	6.9
Child PanoDual	20.6 x 20.5	524 x 521	4.25	1.9	—	—	10.38	4.7	15.12	6.9
Pano Cape*	24.0 x 22.3	610 x 565	2.60	1.2	—	—	4.25	1.9	6.50	3.0
Technician	23.8 x 36.3	603 x 922	5.25	2.4	—	—	7.69	3.5	12.31	5.6
Protective Collar*	5.3 x 27.0	133 x 686	0.50	0.2	—	—	0.94	0.4	1.31	0.6

All sizes and weights are approximate. Length is measured from the neckline, except (*).
Add box weight of .65 lb. (.3 kg) to determine ship weight of aprons [.3 lb (.14 kg) for collars].

Care and Maintenance

Clean with soap and water as needed. Do not use alcohol-base cleaners.
Store at room temperature. Do not store below 50°F (10°C) or above 100°F (37°C).
For storage, hang aprons or lay flat. Do not fold.
Dispose of lead-lined aprons in accordance with federal, state and local laws.



662052 Apron Hanger

The perfect way to store aprons when not in use. Heavy-duty, white epoxy-coated steel. Mounting hardware is included.

*Dimensions: 1.25" high x 14.875" wide x 1.375" deep
(32 x 378 x 35mm)*

DENTSPLY
RINN

A Division of Dentsply International
1212 Abbott Drive
Elgin, Illinois 60123
800 323 0970
847 742 1140
www.rinncorp.com