

Soothe-Guard®

PROTECTIVE RADIOGRAPHY APRONS

Soothe-Guard® Plus⁺

PROTECTIVE RADIOGRAPHY APRONS

Soothe-Guard AIR®

LEAD-FREE RADIOGRAPHY APRONS

S E L E C T I O N G U I D E



DENTSPLY
RINN

Soothe-Guard®

PROTECTIVE RADIOGRAPHY APRONS

The most comfortable aprons we've ever offered!

- Unique silky-smooth material both front and back creates a relaxing and comfortable patient experience
- Smooth draping material with exceptional durability
- Easy to clean — just wipe with damp cloth or disinfectant wipe
- Improved collar design fits better for less gapping, increased patient comfort
- Sturdy sewn-in hanging loops for easy storage
- Velcro closure keeps apron closed securely



.3MM LEAD EQUIVALENCY	Buttercup	Navy	Dove Grey	Light Blue	Tan/Beige	Dusty Mauve
Adult Apron with Collar	668050	669050	6610050	6611050	6612050	6613050
Adult Apron	668049	669049	6610049	6611049	6612049	6613049
Child Apron with Collar	668047	669047	6610047	6611047	6612047	6613047
Child Apron	668048	669048	6610048	6611048	6612048	6613048
Protective Collar	668051	669051	6610051	6611051	6612051	6613051
Pano-Cape Apron	668053	669053	6610053	6611053	6612053	6613053
.5MM LEAD EQUIVALENCY						
Adult Apron with Collar	568050	569050	5610050	5611050	5612050	5613050
Adult Apron	568049	569049	5610049	5611049	5612049	5613049
Child Apron with Collar	568047	569047	5610047	5611047	5612047	5613047
Child Apron	568048	569048	5610048	5611048	5612048	5613048
Protective Collar	568051	569051	5610051	5611051	5612051	5613051
Pano-Cape Apron	568053	569053	5610053	5611053	5612053	5613053

All .5mm lead equivalency aprons are made-to-order (not returnable). 4 to 6 weeks lead time.



Soothe-Guard Protective Collar

Add thyroid protection to any collarless apron.



Soothe-Guard AIR[®]

LEAD-FREE RADIOGRAPHY APRONS



All the features of our regular Soothe-Guard aprons, plus...

- Innovative lead-free alloy lining provides protection equal to lead
- Lining offers up to 25% lighter apron than traditional lead lining
- Easier handling for dental professionals
- Environmentally friendly with no special disposal requirements
- Comfortable Soothe-Guard material
- Improved collar design fits better for less gapping, increased patient comfort
- Sturdy sewn-in hanging loops for easy storage
- Velcro closure keeps apron closed securely



.3MM LEAD EQUIVALENCY	Buttercup	Navy	Dove Grey	Light Blue	Tan/Beige	Dusty Mauve
Adult Apron with Collar	868050	869050	8610050	8611050	8612050	8613050
Adult Apron	868049	869049	8610049	8611049	8612049	8613049
Child Apron with Collar	868047	869047	8610047	8611047	8612047	8613047
Child Apron	868048	869048	8610048	8611048	8612048	8613048
Protective Collar	868051	869051	8610051	8611051	8612051	8613051
Pano-Cape Apron	868053	869053	8610053	8611053	8612053	8613053
Adult Pano-Dual Apron	868054	869054	8610054	8611054	8612054	8613054
Child Pano-Dual Apron	868055	869055	8610055	8611055	8612055	8613055
Technician Apron	868056	869056	8610056	8611056	8612056	8613056

.5MM LEAD EQUIVALENCY

Adult Apron with Collar	878050	879050	
Child Apron with Collar	878047	879047	
Adult Pano-Dual Apron	878054	879054	

All .5mm lead equivalency aprons are made-to-order (not returnable). 4 to 6 weeks lead time.



662052 Apron Hanger

The perfect way to store aprons when not in use. Heavy-duty, white epoxy-coated steel. Mounting hardware is included.

Dimensions: 1.25" high x 14.875" wide x 1.375" deep (32 x 378 x 35mm)

Care and Maintenance

Clean with soap and water as needed. **Do not use alcohol-base cleaners.**

Store at room temperature. Do not store below 50°F (10°C) or above 100°F (37°C).

For storage, hang aprons or lay flat. **Do not fold.**

Dispose of lead-lined aprons in accordance with federal, state and local laws.

New! Soothe-Guard[®] Plus⁺

PROTECTIVE RADIOGRAPHY APRONS

60% more lead in the collar

- Soothe-Guard Plus aprons provide .5mm thyroid protection and .3mm body protection
- Comfortable and easy to clean Soothe-Guard material on both sides
- Collar fastens with Velcro and easily adjusts to patient neck size



.5MM LEAD COLLAR
.3MM LEAD BODY



	Buttercup	Navy	Dove Grey	Light Blue	Tan/Beige	Dusty Mauve
Adult Apron with Collar	468050	469050	4610050	4611050	4612050	4613050
Child Apron with Collar	468047	469047	4610047	4611047	4612047	4613047

All Soothe-Guard Plus aprons have navy blue collars.

Learn more about radiation safety at
www.DentsplyLearning.com/training.asp

Apron Dimensions and Weights

Description	Dimensions		Weight									
	inches	mms	Soothe-Guard Air Lead-Free				Soothe-Guard Lead-Lined				Soothe-Guard Plus	
			.3mm Lead Equiv. lbs.	.3mm Lead Equiv. kgs.	.5mm Lead Equiv. lbs.	.5mm Lead Equiv. kgs.	.3mm Lead Equiv. lbs.	.3mm Lead Equiv. kgs.	.5mm Lead Equiv. lbs.	.5mm Lead Equiv. kgs.	lbs.	kgs.
Adult Patient	23.4 x 25.4	594 x 645	3.8	1.7			4.9	2.2	8.2	3.7		
Adult Patient w/Collar	23.0 x 25.0	584 x 635	4.1	1.8	6.6	3.0	5.5	2.5	8.8	4.0	5.6	2.5
Child Patient	20.9 x 19.7	531 x 500	2.8	1.3			3.3	1.5	5.3	2.4		
Child Patient w/Collar	20.5 x 19.6	521 x 498	2.7	1.2	4.0	1.8	3.6	1.6	5.7	2.6	3.8	1.7
Adult PanoDual	23.3 x 26.8	591 x 679	6.8	3.1	10.4	4.7						
Child PanoDual	20.6 x 20.5	524 x 521	4.3	1.9								
Pano Cape*	24.0 x 11.25	610 x 286	2.6	1.2			3.6	1.6	5.7	2.6		
Technician	23.8 x 36.3	603 x 922	5.3	2.4								
Protective Collar*	5.3 x 27.0	133 x 686	0.5	0.2			0.6	0.3	0.9	0.4		

All sizes and weights are approximate. Length is measured from the neckline, except (*).

Add box weight of .65 lb. (.3 kg) to determine ship weight of aprons [.3 lb (.14 kg) for collars].

DENTSPLY
RINN

A Division of DENTSPLY International Inc.
1212 Abbott Drive
Elgin, Illinois 60123 USA
800-323-0970
847-742-1140
www.RinnCorp.com

